



Prosthodontic

Wax & Modelling Carvers

Cement Spatulas

Wax Knives & Spatulas

Impression Trays

Crown Removers




WAX AND MODELLING CARVERS




Zahle	Zahle	Zahle	Hylin	Hylin	Hylin	Gritman	Lecron	Lecron
TRT-WCZ1	TRT-WCZ2	TRT-WCZ3	TRT-WCH1	TRT-WCH2	TRT-WCH3	TRT-WCG	TRT-WCL1	TRT-WCL4

PKT CARVERS COLOUR CODED



PKT CARVERS



WAX KNIVES & SPATULAS



Fahnenstock
TRT-WKFS



Lessmann
TRT-WKLS



Fahnenstock
TRT-WKFL



Lessmann
TRT-WKLL



Knife
TRT-OK5A



Gritman
TRT-OKG



Flexible
TRT-PS6




Spatula
TRT-PS6



Spatula
TRT-PS4

HEIDEMANN SPATULA



TRT-PFIHS1 #1 #6
TRT-PFIHS2 #1 #6
TRT-PFIHS3 #1 #6

CEMENT SPATULA



TRT-CS6 #3
TRT-CS5 #3
TRT-CS4 #3 #4
TRT-CS3 #3

GINGIVAL CORD PACKER



Serrated
TRT-GCP113
#3 #4

Non-Serrated
TRT-GCP113N
#3 #4

Serrated
TRT-GCPA113
#3 #4

Serrated
TRT-GCPR7
#3 #4

Non-Serrated
TRT-GCPR7N
#3 #4

IMPRESSION TRAYS DENTULOUS PERFORATED

UPPER



TRT-ITRLDPU1

TRT-ITRLDPU2

TRT-ITRLDPU3

TRT-ITRLDPU4

LOWER



TRT-ITRLDPL1

TRT-ITRLDPL2

TRT-ITRLDPL3

TRT-ITRLDPL4

DENTULOUS NON-PERFORATED

UPPER



TRT-ITRLDNPU1

TRT-ITRLDNPU2

TRT-ITRLDNPU3

TRT-ITRLDNPU4

LOWER



TRT-ITRLDNPL1

TRT-ITRLDNPL2

TRT-ITRLDNPL3

TRT-ITRLDNPL4

IMPRESSION TRAYS EDENTULOUS PERFORATED

UPPER



TRT-ITRLEDPU1

TRT-ITRLEDPU2

TRT-ITRLEDPU3

TRT-ITRLEDPU4

LOWER



TRT-ITRLEDPL1

TRT-ITRLEDPL2

TRT-ITRLEDPL3

TRT-ITRLEDPL4

EDENTULOUS NON-PERFORATED

UPPER



TRT-ITRLEDNPU1

TRT-ITRLEDNPU2

TRT-ITRLEDNPU3

TRT-ITRLEDNPU4

LOWER



TRT-ITRLEDNPL1


TRT-ITRLEDNPL2

TRT-ITRLEDNPL3

TRT-ITRLEDNPL4

IMPRESSION TRAY DENTULOUS PERFORATED PEEDO

UPPER



TRT-ITRLDPPU1

TRT-ITRLDPPU2

TRT-ITRLDPPU3

TRT-ITRLDPPU4

LOWER



TRT-ITRLDPPL1

TRT-ITRLDPPL2

TRT-ITRLDPPL3

TRT-ITRLDPPL4

IMPRESSION TRAYS HALF

RIGHT



TRT-ITRLDPHR1

TRT-ITRLDPHR2

TRT-ITRLDPHR3

TRT-ITRLDPHR4

LEFT



TRT-ITRLDPHL1

TRT-ITRLDPHL2

TRT-ITRLDPHL3

TRT-ITRLDPHL4

ANTERIOR AND POSTERIOR



TRT-ITRLDPUA



TRT-ITRLDPLA

ADJUSTABLE



TRT-ITRLDPA



TRT-ITRLDNPA

CALIPER IWANSON



Metal
TRT-CLP1



Wax
TRT-CLP2

ARTICULATING PAPER HOLDER



Miller Straight
TRT-APF2E



Miller Curved
TRT-APF2C



Impression Trays Dentulous Perforated
TRT-ITRLDPP8



Impression Trays Edentulous Perforated
TRT-ITRLEDPP8



Impression Trays Dentulous Non-Perforated
TRT-ITRLDNPP8



Impression Trays Edentulous Non-Perforated
TRT-ITRLEDNPP8

CROWN REMOVER

CROWN REMOVER

TURKEY PATTERN



TRT-CRT TRT-CRT1 TRT-CRT2 TRT-CRT3 TRT-CRT4 TRT-CRT5 TRT-CRT6 TRT-CRT8 TRT-CRT9

MANUAL PATTERN



TRT-CRMP TRT-CRMP1 TRT-CRMP2 TRT-CRMP3

CROWN REMOVER

AUTOMATIC



TRT-CRS



TRT-CRP

Crown Remover Tips with CRS & CRP



TRT-CR1 TRT-CR2 TRT-CR3 TRT-CR4 TRT-CR5 TRT-CR6 TRT-CR7 (for TRT-CRS) TRT-CR8 (for TRT-CRP) TRT-CR9 TRT-CR10

Basic Kit Includes :

- Crown and Bridge Remover # TRT-CRS/TRT-CRP
- Single crown remover tip # TRT-CR1
- Crown remover tip Bifid # TRT-CR2
- Crown & bridge removing tip Sickle # TRT-CR3
- Single crown remover Contra angle tip # TRT-CR4
- Single crown remover tip Inverted curved # TRT-CR5
- Single Crown remover tip straight # TRT-CR6
- Bridge remover Flexi-wire set # TRT-CR7 (for TRT-CRS)
- Bridge remover tip with lock # TR8-CR7 (for TRT-CRP)
- Tips adjustment Key # TRT-CR9
- Open-end adjustment wrench # TRT-CR10

Terms of Use

Removal of :

- Permanently cemented crowns and/or bridges
- Temporarily cemented provisional or final restorations
- Maryland bridges
- Orthodontic brackets/bands
- Removable appliances "stuck" at attachments
- Multi-unit restoration "stuck" on working lab model

Precautions :

1. To remove crowns placed over pins or screw type posts, use lowest "strike intensity" to prevent fracture.
2. This instrument could cause damage to the margins of soft gold crowns. If re-use is anticipated, cut a small groove parallel to the lingual margin.
3. Highest strike intensity could cause root fracture; be careful in long span restorations.
4. Only low intensity strike force should be used with teeth having high mobility and little bone support.
5. For highest strike intensity settings, let the individual situation and common sense prevail.

ARTICULATOR



TRT-A5222

CROWN HOLDING
FORCEP



TRT-CHFERS
12.5cm



TRT-CHFRC
12.5cm

PUTTY CUTTER



PLASTER SHEARS



TRT-PCS
21cm

CROWN OPENER



TRT-CR4566

TEMPORARY
CROWN REMOVING PLIER



TRT-CR4562